


# The handle-needle assembly

**QUICK-RELEASE BUTTON**

disengages the handle to allow reapplication for the second needle passage.

**ERGONOMICALLY-DESIGNED HANDLE & NEEDLE**

are specifically balanced to optimize needle control throughout the entire route of passage.

**CAREFULLY-DESIGNED 4mm-DIAMETER NEEDLES**

offer appropriate rigidity with minimal needle deflection. The 110° arc is designed to follow the trajectory along the back of the pubic bone.

**CONICAL NEEDLE TIP**

is designed to provide smooth needle passage to minimize tissue trauma.

**ALIGNMENT MARKINGS**

provide a visual reference during needle passage.